



\$2,495

SZII MICROSCOPE SPECIALS

MICROSCOPE SOLUTIONS, INC.
 5241 STRATFORD CT.
 CAPE CORAL, FL. 33904
 CELL 949-874-3021/FAX413-460-1119
donwelter@comcast.net
WWW.MICROSCOPESOLUTIONS.NET



call

SZII STEREO ZOOM MICROSCOPE SPECIFICATIONS

Stereo image with longer depth of field & higher resolution

Magnification: 8X~50X

Zoom ratio: 1:6.3

Zoom knob: with detents

Eyepiece inclination: 45°

Interpupillary distance: 52mm~75mm



Diopter: ±5 Diopter

Working distance: 115mm (1X)

Airtight design protects from dust and other contaminants

- Antistatic:

5 Year Warranty

		<u>PART NO</u>	<u>DESCRIPTION</u>	<u>PRICE</u>
1	1	SZ2-200	SZ-II Zoom Stereo Binocular Microscope with a Track Stand, including SZ401101 Stereo Zoom Binocular Body, SZ402101 10X Adjustable Eyepieces (Pair), MA221101 76mm Track Stand, SZ403201 49.5mm Ring Adapter	\$795.00
				
2	1	SZ2-1000	SZ2 Zoom Stereo Binocular Microscope with a Boom Stand, including SZ401101 Stereo Zoom Binocular Body(.8-5X), SZ402101 10X Adjustable Eyepieces (Pair), MA231101 Boom Stand, & MA251101 Focus-Arm.	\$1,095.00
				
3	1	MA312101	Fiber Optic Illuminator,21V 150W Halogen light bulb(EKE), 5/8"	\$ 220.00
			Fiber Optic Cable adapter with metal housing.	
4	1	MA322101	Dual Pipe Light Guide with Focus Lens	\$ 150.00
5	1	MA321101	Annular Ring Light Guide, Ring's I.D.: 58mm; 1meter in length	\$ 270.00
6	7	SZ403201	SZ403201 49.5mm Ring Adapter	\$ 12.00

